# 3D Accuitomo 170

# 3D Accuitomo 170

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
The i-Dixel imaging software offers a wide variety of features to help you quickly and easily create comprehensive treatment plans and explain those plans to your patients. Mandibular canal marking, implant presentation, multiplanar reconstruction are just a few of the features that i-Dixel provides for diagnoses. i-Dixel is also fully DICOM compliant and provides quick and easy integration with both practice management software and advanced treatment planning tools.

## Volume Rendering


Volume rendering displays a solid 3D image showing the bone structure and dental arch. The volume rendered image is linked to the slices and moves in real time whenever a slice is adjusted.

From multiplanar reconstructions to implant presentations

## Implant Presentation

Select from a variety of popular implant manufacturers and place the implant directly in the CT slice for presentation. This presentation can be easily understood by patients and helps with case acceptance.


# i-Dixel WEB


## Freedom from platform and simplicity of design

i-Dixel WEB runs as a web service on an X-ray server PC included with your Morita X-ray system. It serves as a local and secure web-based dental image processing service that you can access throughout your practice on a wide range of devices.

## No software installation needed

With the latest advancements in web technology, i-Dixel WEB gives you the freedom to view your images wherever you want and however you want. Gone are the days of complicated chairside PC setups, and limited choices of hardware. Mac OS X and even iPads can be used to view and edit data from a Morita X-ray system.

## 3D image processing

Even though i-Dixel WEB is a fully web-based system, the features you've come to expect are not compromised. View Morita's signature high quality images on the device of your choice!

Mac, iPad, Safari, OS X and Mac OS are trademarks of Apple Inc., registered in the U.S. and other countries. Windows, Windows 8, Windows 10 and Internet Explorer are registered trademarks of Microsoft Corporation in the United States and other countries.

# Implantology

Image 1a: Clinical aspect at the initial examination

Case 1: Female patient referred for 3-dimensional analysis of an esthetic complication in the left maxillary incisor region (a).

The clinical status exhibits a mucosal recession as well as a flattening and discoloration of the facial mucosa at the implant crown. The patient complained about recurrent peri-implant infections.

\*: Nasal palate tube

Image 1b: Coronal CBCT slice

Image 1c: Axial CBCT slice

Image 1d: Sagittal CBCT slice

Case courtesy Prof. em. Dr. Daniel Buser, Clinic for oral surgery and somatology University Bern Prof. Dr. Michael Bornstein, Department of Oral Health & Medicine, University Center of Dental Medicine Basel UZB, University of Basel (Switzerland) and Honorary Professor in Oral and Maxillofacial Radiology, Faculty of Dentistry, The University of Hong Kong (Hong Kong SAR, China).

# Periodontics

Image 2a: Para-coronal view


Image 2b: Cross-sectional view


Image 2c: Axial view


Image 2d: Volumetric Rendered view

An asymptomatic 51-year-old male presents with numerous complaints associated with the mandibular anterior teeth including unesthetic “black triangles” between the teeth, shrinking gums making the teeth look long, and “loose”. The patient reports an extensive dental history of periodontal bone loss and therapy. A 4 cm x 4 cm FOV at 0.08 nominal voxel size was acquired and para-coronal (a), cross-sectional (b), axial (c), and three dimensional volumetric rendered (d) images clearly show the presence of calculus on the mandibular central incisors and generalized alveolar peri-circumferential radiographic bone loss extending to the middle third of the root and beyond consistent with a diagnosis of Periodontitis Stage III, localized incisor pattern. The patient was referred to a periodontist for management.

Case courtesy, Drs William C. Scarfe and Gustavo Santaella, Louisville Kentucky

# Endodontics


Image 3a: Cross-sectional view

Image 3b: Para-sagittal view

An asymptomatic 27-year-old female presents immediately after fixed orthodontic appliance therapy with bitewing imaging showing diffuse radiolucent increase in the middle third of the pulp chamber of the maxillary right canine compared to the contralateral side. An 8 cm x 8 cm FOV at 0.25 nominal voxel size was acquired and a suspicious opacification of the pulp canal of the maxillary right canine noted on axial (a) images. A 4 cm diameter "zoom reconstruction" was performed at a nominal 0.08 mm voxel resolution centered on the maxillary canine. Cross-sectional (a), and para-sagittal (b) images clearly identify dentin resorption on the mesio-palatal aspect of the root extending peri-circumferentially around the pulp canal and chamber and loss of the root continuity at the level of the cervical third of the root indicative of Class III invasive cervical resorption. The patient was referred to an endodontist for management.

Image 4a: Axial view

Image 4b: Cross-sectional view

A 24-year-old female presents with vague tooth sensitivity in the maxillary right region. Her dental history reveals active orthodontics more than 8 years previously involving extrusion of her right maxillary canine. Intraoral imaging is noncontributory. An 8 cm x 8 cm FOV at 0.25 nominal voxel size was acquired and a suspicious opacification of the pulp canal of the maxillary right canine noted. A 4 cm diameter "zoom reconstruction" was performed at a nominal 0.08 mm voxel resolution centered on the maxillary canine and these indications are noted on axial (a) and magnified cross-sectional (b) images. This image serves as a baseline for periodic, limited field, high resolution CBCT imaging to determine progression.

# Oral Surgery

Image 5a: MPR Panoramic view

Image 5b: Para-sagittal view

An asymptomatic 56 year-year-old male presents with a history of incidental discovery of possible mandibular pathology in the left mandible on routine panoramic imaging taken 1 month previously. A 10 cm x 10 cm FOV at 0.25 nominal voxel size was acquired. A reformatted MPR panoramic (a), and para-sagittal (b) images clearly show a single, well-defined, corticated, irregularly-shaped bilobular low density lesion anterior to the lingula and mandibular foramen within the left ascending ramus with extension inferiorly through the intramedullary bone to include the entire alveolus posterior to the distal root of the left mandibular third molar. There is lingual cortical expansion but no perforation. The mandibular canal is intact throughout its course within the ramus and the mandible. The patient was referred to an oral surgeon and excisional biopsy reveals an odontogenic keratocyst.

Case courtesy, Drs William C. Scarfe and Gustavo Santaella, Louisville Kentucky

Image 6a: Axial view

Image 6b: Volume Rendered view

Buccal bifurcation cyst. Buccal localization of the lesion with an extended resorption of the vestibular cortical plate and a periosteal reaction.


Case courtesy, Prof. Dr. Reinhilde Jacobs, Center for Dentomaxillofacial Imaging, University Hospitals Leuven and Department of Imaging & Pathology, OMFS-IMPATh Research Group, KU Leuven, Leuven, Belgium

# Specifications

Trade Name:	3D Accuitomo XYZ Slice View Tomograph
Model:	MCT-1
Type:	EX1/2 F17
Power Supply:	AC 100/ 110/ 120 V AC 220/230/240 VAC
Power Consumption:	max 2.0 kVA
Dimensions:	
Main Unit:	W1,620 mm x D1,250 mm x H2,080 mm (63-3/4" x 49-1/4" x 82")
Control Box:	W100 mm x D40 mm x H115 mm (4" x 1-5/8" x 4-1/2")
Weight:	Approx. 400kg (Approx. 882lbs)
X-ray Head	
Tube Voltage:	60-90 kV
Tube Current:	1-10 mA (Max 8mA : Hi-Fi, Hi-Res Mode)
Focal Spot Size:	0.5
Exposure Time:	Std Mode: 17.5 / 9.0 sec Hi-Fi Mode : 30.8 / 15.8 sec Hi-Res Mode : 30.8 / 15.8 sec Hi-Speed Mode : 10.5 / 5.4 sec
Field of View	ø 40 x H40 mm ø 60 x H60 mm ø 80 x H80 mm ø 100 x H50 mm ø 100 x H100 mm ø 140 x H50 mm ø 140 x H100 mm ø 170 x H120 mm
Voxel Size	80µm / 125µm / 160µm / 250µm

\* X-ray protection should be provided for the patient when X-rays are emitted.  
\* Design and specifications are subject to change without notification.

Development and Manufacturing

**J. MORITA MFG. CORP.**

680 Higashihama Minami-cho, Fushimi-ku,  
Kyoto 612-8533, Japan  
T +81. (0)75. 611 2141, F +81. (0)75. 622 4595

**Morita Global Website**

**[www.morita.com](http://www.morita.com)**

Distribution

**J. MORITA CORP.**

3-33-18 Tarumi-cho, Suita-shi, Osaka 564-8650, Japan  
T +81. (0)6. 6380 1521, F +81. (0)6. 6380 0585

**J. MORITA USA, INC.**

9 Mason, Irvine CA 92618, USA  
T +1. 949. 581 9600, F +1. 949. 581 8811

**J. MORITA EUROPE GMBH**

Justus-von-Liebig-Strasse 27b, 63128 Dietzenbach, Germany  
T +49. (0)6074. 836 0, F +49. (0)6074. 836 299

**MORITA DENTAL ASIA PTE. LTD.**

150 Kampong Ampat  
#06-01A KA Centre, Singapore 368324  
T +65. 6779. 4795, F +65. 6777. 2279

**J. MORITA CORP. AUSTRALIA & NEW ZEALAND**

Suite 2.05, 247 Coward Street, Mascot NSW 2020, Australia  
T +61. (0)2. 9667 3555, F +61. (0)2. 9667 3577

**J. MORITA CORP. MIDDLE EAST**

4 Tag Al Roasaa, Apartment 902, Saba Pacha 21311 Alexandria, Egypt  
T +20. (0)3. 58 222 94, F +20. (0)3. 58 222 96

**J. MORITA CORP. INDIA**

Filix Office No.908, L.B.S. Marg, Opp. Asian Paints, Bhandup (West), Mumbai 400078, India  
T +91-22-2595-3482

**J. MORITA MFG. CORP. INDONESIA**

28F, DBS Bank Tower, Jl. Prof. Dr. Satrio Kav. 3-5, Jakarta 12940, Indonesia  
T +62-21-2988-8332, F + 62-21-2988-8201

**SIAMDENT CO., LTD.**

71/10 Mu 5, Thakham, Bangpakong, Chachuengsao 24130, Thailand  
T +66. 38. 573042, F +66. 38. 573043  
[www.siamdent.com](http://www.siamdent.com)

Diagnostic and Imaging Equipment

Treatment Units

Handpieces and Instruments

Endodontic System

Laser Equipment

Laboratory Devices

Educational and Training Systems

Auxiliaries

Subject to technical changes and errors.